# Navodila za uporabo

## Instructions for use

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Thank you for choosing our product. The NOEL Žebrowski sp. j. has designed and manufactured Base Stabilo System to the highest safety and functional use standards. We are confident that our product will meet your highest requirements and expectations. These instructions for use are an integral part of vacuum positioning system stabilizing torso, limbs and head Base Stabilo System. Please read the instructions thoroughly before using the product.

#### Sistem za vakuumsko pozicioniranje, ki stabilizira trup, okončine in glavo Base Stabilo System

#### Vacuum positioning system stabilizing torso, limbs and head Base Stabilo System

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
Vacuum Positioning Base Stabilo System is a product designed for children and adults with cerebral palsy and myelomeningocele, compensating for and alleviating the disability by stabilising hips, pelvis and thighs in sitting position. Positioning Base Stabilo System can be used for alleviating the effects of remaining in one position for extended periods (e.g. contractures).

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
Base Velcro Stabilo System cushions are covered a special outer material that allows attaching velcro strips. Base Hilo Velcro Stabilo System has an additional abduction wedge that reduces contractures and ensures proper placement of the legs in a sitting position. Base Doro has side laterals which tuck whole pelvis up, and the neoprene strap to increase stability. Doro enables hips correction, stabilization and relieving its pressure at the same time.

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#### PRODUCT SAFETY WARNINGS

- positioning Base Stabilo System must be securely fixed to seats or backrests
- use of airy cover is advisable if increased ventilation is required around the body
- avoid contact with sharp objects that may damage the cushion
- do not clean using solvent-based cleaners, e.g. toluene or acetone do not pull on valve
- do not overinflate cushion
- use only special pump included with cushion for inflation
- do not iron
- in case of tear, micro-granules may escape which may be dangerous for small children
- only persons who received instructions or training may use the device

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
Permitted user weight 80 kg

**ATTENTION**

- the thickness of the cushion can naturally reduce by 5% in the initial period of use
- in case of considerable changes of temperature use pump again to obtain proper hardness (higher temperature softens the cushion, lower temperature hardens the cushion)

**INSTRUCTIONS FOR USE**

Hand pump

1. Take out tube in pump handle.
2. Screw tube into cushion valve and wait until cushion reaches the required softness.
3. Shape cushion (either manually to adjust shape or by sitting or leaning on it).
4. Screw pump into tube and deflate air until desired hardness has been reached.
5. Remove tube from pump and screw protecting cap on valve. To
6. change shape - screw tube into valve and wait until cushion softens, then repeat steps from 2 to 5.

Electric pump

1. Take out tube from box with electric pump.
2. Screw off valve cap, screw tube (tapped part) into cushion valve, and push second tube into SOFT socket in electric pump. Turn on pump and wait until cushion will be soft enough, then turn pump off.  
Start shaping cushion (either manually to adjust shape or by sitting or leaning on it).
3. Change position of tube from SOFT to HARD and turn pump on
4. removing air until the desired hardness has been reached. Screw off tube, put cap back on and tighten.
5. To change the shape repeat steps from 2 to 5.

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**MAINTENANCE**

- regularly check condition of product – pay attention to any cracks and rips
- screw on valve cap before washing
- hand wash cushion under running water using a standard detergent, e.g. a washing up liquid
- suitable for washing in washing machine at 40°C cycle. Fabric softener is sufficient, does not require use of washing powder wipe dry with paper towel
- standard disinfectants used in hospitals can be used
- do not use solvent-based cleaners such as toluene or acetone
- for Velcro products, Hoover the soft surface regularly to prevent loss of adherence

**WARRANTY**

- the warranty is provided for a period of 24 months
- the warranty does not cover mechanical damage to Base Stabilo System
- if in doubt, please contact seller or manufacturer directly report any incidents to seller or manufacturer directly

**DECLARATION OF CONFORMITY**

Vacuum positioning system stabilizing torso, limbs and head Base Stabilo System medical device complies with the basic requirements of the Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

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